

RB-PEG-SH (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	RB-PEG-SH (MW 2000) (Rhodamine B-PEG-SH (MW 2000)) is a fluorescent dye comprised of Rhodamine B, PEG, and a thiol (SH) group. The highly reactive SH can form stable thioether bonds with various functional groups, such as maleimide, exhibiting excitation/emission spectra of 546/610 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481