

P2Y12R ligand-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	P2Y12R ligand-1 (Compound 41) is a high-affinity P2Y12 receptor ligand with a K_i of 17.1 nM. The [18F]-labeled version of P2Y12R ligand-1 is suitable for research involving brain PET imaging.
Targets(IC50)	P2Y Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481