

Pam3CSK4-Biotin TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pam3CSK4-Biotin (Pam3Cys-Ser-(Lys)4-Biotin) TFA is the biotinylated version of Pam3CSK4. It acts as an agonist for Toll-like receptor 1/2 (TLR1/2).
Targets(IC50)	TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481