

4Arm-PEG-RB (MW 20000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | 4Arm-PEG-RB (MW 20000) [4Arm-PEG-Rhodamine B (MW 20000)] is a fluorescent dye with a four-arm structure made of Rhodamine B and PEG. It is suitable for applications in drug delivery and fluorescent imaging (Ex/Em = 546/610 nm). |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481