

## D-myo-Inositol-2,3,5,6-tetraphosphate tetrasodium

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	D-myo-Inositol-2,3,5,6-tetraphosphate (2,3,5,6-IP4) tetrasodium is an inositol phosphate utilized in the study of calcium channels.
Targets(IC50)	Calcium Channel

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481