

RB-PEG-COOH (MW 3400)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RB-PEG-COOH (MW 3400) (Rhodamine B-PEG-COOH (MW 3400)) is a fluorescent dye comprising Rhodamine B, PEG, and COOH. It is suitable for fluorescence labeling and imaging, with excitation and emission wavelengths of 546 nm and 610 nm, respectively.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481