

RB-PEG-Biotin (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	RB-PEG-Biotin (MW 10000) (Rhodamine B-PEG-Biotin (MW 10000)) is a fluorescent dye composed of Rhodamine B, PEG, and biotin. Biotin binds with high affinity to streptavidin or avidin. RB-PEG-Biotin (MW 10000) plays a crucial role in the specific capture and detection of biomolecules (Ex/Em = 546/610 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481