

GSNOR-IN-2

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | GSNOR-IN-2 (Compound 45) is an inhibitor of S-nitrosoglutathione reductase (GSNOR) with an IC50 of 44.12 nM. It exhibits neuroprotective properties and shows potential for research in ischemic stroke. |
| Targets(IC50) | GSNOR  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481