

m-PEG-NHS ester (MW 30000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | m-PEG-NHS ester (mPEG-SC) (MW 30000) is useful for modifying active molecules to enhance their antigenicity and immunogenicity, aiding in the formulation of injectable preparations. Modification of serine protease lumbrokinase (LK) with m-PEG-NHS ester does not affect its fibrinolytic and thrombolytic activities, indicating promising application potential. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481