

(D-Trp32)-Neuropeptide Y (porcine)

Chemical Properties

CAS No. : 153549-78-1

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(D-Trp32)-Neuropeptide Y (porcine) is a highly selective agonist of the neuropeptide Y (NPY) Y5 receptor. It exhibits appetite-stimulating activity and can inhibit Forskolin-induced cAMP formation.
Targets(IC50)	Neuropeptide Y Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481