

5-cis-15(R)-Iloprost

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	5-cis-15(R)-Iloprost (5-cis-15(R)-Ciloprost) is the C-5 cis-isomer and 15(R)-epimer isomer of Iloprost. Iloprost, a prostacyclin (PGI ₂) analog, plays a role in embryonic development, alleviates inflammation, and inhibits tumor metastasis. Additionally, Iloprost is applicable in peripheral vascular research.
Targets(IC50)	Others,Prostaglandin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481