

CA-II/Dkk1-IN-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CA-II/Dkk1-IN-1 (Compound 5d) is a dual inhibitor of CA-II and Dkk1, with an IC50 of 6.90 nM for CA-II. It exhibits significant antioxidant activity and exceptional DNA binding capability. CA-II/Dkk1-IN-1 is applicable in cancer research, including esophageal cancer, renal cancer, and lung cancer.
Targets(IC50)	Carbonic Anhydrase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481