

FA-PEG-COOH (MW 2000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FA-PEG-COOH (MW 2000) (Folate-PEG-COOH (MW 2000)) is a PEG derivative modified with folic acid, featuring a carboxyl group (-COOH). The carboxyl group readily forms stable amide bonds with amino groups and ester bonds with hydroxyl groups. Folic acid (-FA) has a high affinity for folate receptors, making it suitable for cell membrane receptor-targeted drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481