

4-Arm-PEG-DA (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4-Arm-PEG-DA (MW 2000) [also known as 4-Arm-PEG-Acrylate (MW 2000)] is a polymer featuring a four-arm PEG core with acrylate groups chemically bonded to the end of each PEG chain. This structure is suitable for drug delivery applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481