

RB-PEG-DSPE (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	RB-PEG-DSPE (MW 1000) (Rhodamine B-PEG-DSPE (MW 1000)) is a PEG-phospholipid labeled with Rhodamine B. It is suitable for the preparation of liposomes or nanoparticles intended for the delivery of nutrients such as mRNA or DNA vaccines (Ex/Em = 546/610 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481