

RB-PEG-RB (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | RB-PEG-RB (MW 3400) (Rhodamine B-PEG-Rhodamine B (MW 3400)) is a fluorescent dye composed of two Rhodamine B molecules and polyethylene glycol (PEG). It is suitable for fluorescence labeling and imaging with excitation/emission wavelengths of 546/610 nm. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481