

DNP-PEG3-NHS ester

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DNP-PEG3-NHS ester is a PEG derivative that includes DNP (2,4-dinitrophenyl), three PEG units, and an NHS ester. The NHS ester can form covalent bonds with amino-containing molecules (such as proteins, antibodies, etc.). This compound is applicable in drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481