

6-Arm-PEG-FA (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	6-Arm-PEG-FA (MW 3400) [6-Arm-PEG-Folate (MW 3400)] is a PEG derivative modified with folic acid. Polyethylene glycol enhances solubility and stability, reduces the immunogenicity of peptides and proteins, and inhibits nonspecific binding of charged molecules to modified surfaces. Folic acid (-FA) targets cell membrane receptors for drug delivery due to its high affinity for folate receptors.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481