

Maleimide-Ph(3,5-F)-PEG4-Val-Ala

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Maleimide-Ph(3,5-F)-PEG4-Val-Ala is a cleavable ADC linker belonging to the PEG class, utilized in the synthesis of antibody-drug conjugates (ADCs). |
| Targets(IC50) | ADC Linker |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481