

FMOC-NH-PEG-SH (MW 5000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FMOC-NH-PEG-SH (MW 5000) (FMOC-NH-PEG-Thiol (MW 5000)) is a PEG derivative containing an Fmoc protecting group and a reactive thiol (-SH) functional group. The thiol group is highly active and can specifically react with various molecules to form stable covalent bonds. Fmoc serves as a common amino-protecting group, shielding the amine from unintended reactions until it is intentionally removed for specific coupling reactions.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481