

Bis-PEG-COOH (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bis-PEG-COOH (MW 1000) is a chemically identical functional group disubstituted PEG derivative. The carboxylic acid groups in Bis-PEG-COOH (MW 1000) can react with compounds containing amino or hydroxyl groups to form stable chemical bonds, making it suitable for drug delivery applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481