

YVAD-CHO TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	YVAD-CHO TFA is an inhibitor of interleukin-1 β converting enzyme (ICE, caspase 1) with a K_i of 0.76 nM. It suppresses the production of mature IL-1 β and can partially slow the death of motor neurons in mice with facial nerve injuries.
Targets(IC50)	Caspase,Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481