

Ac-RYQK(Ac)-AMC

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ac-RYQK(Ac)-AMC (Compound S5) is a fluorescent substrate based on the H3 sequence used for studying SIRT6. It mimics the deacetylation site of H3K56, significantly enhancing deacetylation signals while exhibiting an excellent signal-to-noise ratio. This compound is applicable in aging and cancer research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481