

Cipepofol glucuronide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cipepofolglucuronide (Ciprofol Glucuronide) (M4) is a glucuronide metabolite of Ciprofol. This compound does not significantly inhibit P450s or UDP-glucuronosyltransferase (UGT) enzymes and poses a low risk as a perpetrator of drug-drug interactions (DDIs). Cipepofolglucuronide is used in anesthesia research.
Targets(IC50)	Drug Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481