

## Maleimide-PEG-COOH (MW 5000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DMG-PEG-Mal (MW 5000) is a PEG derivative composed of N,N-Dimethylglycine (DMG), PEG units, and maleimide. The maleimide component reacts with thiol groups (-SH) to form stable thioether bonds. DMG-PEG-Mal (MW 1000) can be utilized for drug delivery purposes.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481