

Carbaprostacyclin-biotin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Carbaprostacyclin-biotin (cPGI-biotin; Carbacyclin-biotin) is a Carbacyclin linked to biotin. Carbacyclin is an analog of PGI ₂ , functioning as a prostacyclin (PGI ₂) receptor agonist and a vasodilator, effectively inhibiting platelet aggregation.
Targets(IC50)	Prostaglandin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481