

## PD-L1 inhibitory peptide TFA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PD-L1 inhibitory peptide TFA is an inhibitory peptide targeting programmed cell death ligand 1 (PD-L1). It binds to PD-L1, lifting immune suppression and restoring the antitumor activity of T cells. PD-L1 inhibitory peptide TFA shows potential for use in cancer research.
Targets(IC50)	PD-1/PD-L1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481