

## Azido-PEG3-DYKDDDDK TFA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Azido-PEG3-DYKDDDDK (azide-PEG3-FLAG) TFA serves as a multifunctional fusion tag for the purification of recombinant proteins. This click chemistry reagent comprises an azide group, enabling it to undergo copper-catalyzed azide-alkyne cycloaddition (CuAAC) with alkyne-containing molecules. Additionally, it can participate in strain-promoted alkyne-azide cycloaddition (SPAAC) with molecules that possess DBCO or BCN groups.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481