

6-Arm-PEG-Silane (MW 600)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	6-Arm-PEG-Silane (MW 600) is a silane-modified PEG derivative that can modify surfaces such as glass, silicon, or others through reactions between hydroxyl and ethoxy (methoxy) silane groups. PEG effectively inhibits non-specific surface binding.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481