

RB-PEG-DSPE (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | RB-PEG-DSPE (MW 10000) (Rhodamine B-PEG-DSPE (MW 10000)) is a Rhodamine B-labeled PEG phospholipid. It is utilized in the preparation of liposomes or nanoparticles to deliver nutrients, such as mRNA or DNA vaccines, with excitation/emission wavelengths of 546/610 nm. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481