

RB-PEG-FA (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	RB-PEG-FA (MW 2000) (Rhodamine B-PEG-FA (MW 2000)) is a fluorescent dye composed of Rhodamine B, PEG, and Folic Acid (FA). It is suitable for applications in cell imaging, folate receptor targeting, and detection, with excitation and emission wavelengths of 546/610 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481