

RB-PEG-FA (MW 2000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RB-PEG-FA (MW 2000) (Rhodamine B-PEG-FA (MW 2000)) is a fluorescent dye composed of Rhodamine B, PEG, and Folic Acid (FA). It is suitable for applications in cell imaging, folate receptor targeting, and detection, with excitation and emission wavelengths of 546/610 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481