

## 6-Arm-PEG-Mal (MW 3400)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	6-Arm-PEG-Mal (MW 3400) is a Maleimide-modified PEG derivative that forms stable thioether bonds with sulfhydryl (-SH) groups. The PEG modification enhances solubility and stability and reduces the immunogenicity of proteins and peptides.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481