

DBCO-PEG-LA (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | DBCO-PEG-LA (MW 3400) is a PEG derivative composed of DBCO, PEG units, and Lipoic Acid. This compound contains a DBCO group, enabling it to participate in strain-promoted alkyne-azide cycloaddition (SPAAC) reactions with azide-containing molecules. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481