

## USP7 Ligand-Linker Conjugates 1

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	USP7Ligand-Linker Conjugates 1 is a Target Protein Ligand-Linker Conjugate comprising a USP7 ligand and a PROTAC linker, capable of recruiting E3 ligase. It is utilized in the synthesis of PROTACUSP7 Degradar-1.
Targets(IC50)	DUB,Target Protein Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481