

Exendin-4, Cy5 labeled TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Exendin-4, Cy5 labeled TFA is a Cy5-tagged version of Exendin 4. It serves as an in vivo β -cell tracer. Exendin 4 is a long-acting GLP-1 receptor agonist.
Targets(IC50)	Others,Glucagon Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481