

## 6-Arm-PEG-Silane (MW 2000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

|               |  |
|---------------|--|
| Description   | 6-Arm-PEG-Silane (MW 2000) is a silane-modified PEG derivative designed to modify surfaces such as glass, silicon, or other materials by reacting through hydroxyl groups with ethoxy (or methoxy) silane. PEG functions to inhibit non-specific binding on the surface. |
| Targets(IC50) | Others   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481