

PROTAC BRD4 Degradator-37

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	PROTACBRD4 Degradator-37 is a PROTACBRD4 degrader with a DC50 of 36.4 nM and a Dmax of 73% in PANC-1 cells. It exhibits cytotoxicity against PANC-1 cells with a GI50 of 0.282 μ M. PROTACBRD4 Degradator-37 is applicable for cancer research.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481