

EGFR-IN-47-C6-COOME

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	EGFR-IN-47-C6-COOME is a synthetic target protein ligand-linker complex used in the synthesis of PROTACs, such as PROTACEGFRdegrader 14. This compound, PROTACEGFRdegrader 14, functions as an EGFR PROTAC degrader with antitumor activity.
Targets(IC50)	Target Protein Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481