

Ac-DMQD-AMC TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ac-DMQD-AMC TFA is a fluorescent substrate for Caspase-3 or Caspase-10, analogous to Ac-DEVD-AMC, capable of releasing the fluorescent group AMC (7-amino-4-methylcoumarin), thereby emitting fluorescence. Ac-DMQD-AMC TFA is useful for measuring Caspase activity in vitro during apoptosis monitoring.
Targets(IC50)	Caspase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481