

RB-PEG-RB (MW 6000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RB-PEG-RB (MW 6000) (Rhodamine B-PEG-Rhodamine B (MW 6000)) is a fluorescent dye made up of two Rhodamine B molecules and PEG. It is suitable for fluorescence labeling and imaging with excitation/emission wavelengths of 546/610 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481