

DBCO-PEG16-NHS ester

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DBCO-PEG16-NHS ester is a PEG derivative composed of DBCO, PEG units, and an NHS ester. The NHS ester can react with amino acids or other amine-containing molecules. The DBCO group in DBCO-PEG16-NHS ester can engage in a strain-promoted alkyne-azide cycloaddition reaction (SPAAC) with azide-containing molecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481