

BOC-NH-PEG-FA (MW 3400)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BOC-NH-PEG-FA (MW 3400) [BOC-NH-PEG-Folate (MW 3400)] is a PEG derivative modified with folic acid. The folic acid component (-FA) has a strong affinity for folate receptors, making it suitable for targeting cell membrane receptors in drug delivery
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481