

RGD-DA6

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RGD-DA6 is an $\alpha\beta3$ -targeting ligand composed of RGDfK (an RGD peptide targeting integrin $\alpha\beta3$ ), DOTADG, d-Glu, Lys, AMBA, and IPBA. Through radiolabeling, RGD-DA6 forms an $\alpha\beta3$ -targeted radioligand for high-contrast tumor imaging and improved anticancer efficacy.
Targets(IC50)	Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481