

## 6-Arm-PEG-FA (MW 600)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	6-Arm-PEG-FA (MW 600) (6-Arm-PEG-Folate (MW 600)) is a PEG derivative modified with folic acid. Polyethylene glycol enhances solubility and stability, decreases the immunogenicity of peptides and proteins, and inhibits nonspecific binding of charged molecules to modified surfaces. Folic acid (-FA) has a high affinity for folate receptors, making it suitable for targeted drug delivery to cell membrane receptors.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481