

DBCO-PEG-SH (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DBCO-PEG-SH (MW 2000) (DBCO-PEG-Thiol (MW 2000)) is a PEG derivative modified with a DBCO group. It contains the DBCO moiety, enabling it to undergo a strain-promoted alkyne-azide cycloaddition reaction (SPAAC) with molecules that have azide groups.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481