

## Endo-BCN-PEG5-amine

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Endo-BCN-PEG5-amine (Endo-BCN-PEG5-NH <sub>2</sub> ) is a click chemistry crosslinker. The BCN group exhibits high reactivity with azide-labeled molecules. The amine group reacts effectively with carboxylic acids, activated NHS esters, and carbonyl compounds (ketones, aldehydes), among others.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481