

FMOC-NH-PEG-SH (MW 1000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FMOC-NH-PEG-SH (MW 1000) (FMOC-NH-PEG-Thiol (MW 1000)) is a PEG derivative featuring an Fmoc protective group and a thiol (-SH) reactive group. The thiol group is highly reactive, enabling it to selectively bond with various molecules to form stable covalent links. Fmoc is a common amino protecting group that safeguards against unwanted reactions until a specific coupling reaction is required.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481