

BODIPY 480/508-cholesteryl oleate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	BODIPY 480/508-cholesteryl oleate is a fluorescently labeled derivative of Cholesteryl oleate. It is utilized for monitoring lipid transport mechanisms in mammalian cells (Ex/Em = 480/508 nm).
Targets(IC50)	Others,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481