

RB-PEG-DSPE (MW 2000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RB-PEG-DSPE (MW 2000) (Rhodamine B-PEG-DSPE (MW 2000)) is a PEG phospholipid-labeled Rhodamine B. It can be utilized in the preparation of liposomes or nanoparticles for the delivery of nutrients, such as mRNA or DNA vaccines (Ex/Em = 546/610 nm).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481